



BBR-B Series Double Door

| Blood Bank

☰ Details Info

Specifications:

Model	BBR-650B	BBR-1000B
Capacity	650L	1000L
Temp. Range	4 ±1°C	
Temp. Accuracy	0.1°C	
Control System	Microprocessor control	
Display	Large screen LED display	
Alarm	Audio and visual alarm for: High or low temperature, Sensor failure, System failure; Door ajar; Power failure,	

	Condenser abnormal	
Refrigeration Type	Forced air refrigeration system	
Refrigerant	R134a, CFC-Free	
Defrost	Auto defrost	
Construction	Unibody Design, Mono-assembly foaming (Rigid polyurethane insulation material)	
Internal Material	Stainless steel	
External Material	Cold-rolled steel with anti-bacteria powder coating	
Door	Three-ply transparent glass door with lock, Heating film preventing frosting, Sealing joint strip with magnetic force	
Drawers/Shelves	Adjustable shelves	SS drawers
Casters	4 casters	4 casters
Consumption	364W	600W
Power Supply	AC220V±10%, 50/60Hz	AC220V±10%, 50/60Hz
External Size (L*W*H), mm	1220*648*1885	1260*700*1980
Standard Accessory	Temperature printer, 12 shelves, 36 baskets, fluorescent light	Temperature printer, 12 drawers, fluorescent lamp
Blood Bags(450ml)	360 pcs	576 pcs

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office
: warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved