



## BBT-10 Blood Bag Tube Sealer

| Blood Bank

☰ Details Info

### Features:

Blood Bag Tube Sealer Model BBT-10 , Bag Sealer, Tube Sealing Machine & Hospital Laboratory equipment to seal the blood bag tube By high frequency dielectric heating system. The sealing can be achieved in less then 1.5 seconds per tube by simple operation. Complete sealing of the blood bag tube is annexure by high frequency dielectric heating. Sealing of each tube takes approx 1.5 seconds and is thus very efficient. Pulling the sealing section to both sides can readily cut sealed tube

### Specification:

Oscillation Frequency	40.68 MHZ
Radio Frequency output	15W
Sealing time	Approx 1.5 sec

Input Volt	230V AC
Power Frequency	50Hz
Power Consumption	100 V A
Weight	5.5 Kg approx
Size	180mmx215mmx245mm (HxDxW)



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved