



## BCM - Series Blood Collection Monitor

| Blood Bank

Details Info

Blood Collection Monitor is a compact unit that provide safe and smooth action to ensure Mixing of Blood and Anticoagulant without damaging Blood Cells during collection of Blood from a donor.

### Features:

- Pre-selection of volume Auto tare facility Display of Set Volume
- Process Volume, weight and Time Volume can be set 1 ml
- Increment with two set points Micro control based program Volume can be set from 1 ml. To 999ml
- Provision for pausing collection and change
- Program during blood collection
- Auto Calibration Facility Battery power up to 7 to 8 Hrs. (Only in Battery backup model)

### Specification:

Catalog #	BCM-10	BCM-20
Input Voltage	230 V AC	230 V AC (SMPS)

Oscillation Speed	10 RPM +/- 1	
Display	LCD 16*2	
Alarm	Audio Visual	
Battery Backup	-	7-8 Hrs



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved