



BDC-10 Blood Donor Chair

| Blood Bank

☰ Details Info


Micro-controller and Dual actuator based blood donor couch, Blood Donor Chair, Hospital Chair & Hospital Furniture is specially designed to make blood withdrawals easier, safe and provide comfortable position for the donor. Donor's position can be easily adjusted to form either a comfortable easy chair or a bed. All surfaces are clean and smooth. Easily accessible and can be kept clean. Electrical and mechanical accessories are housed inside the structure making the product safer. Seat elements are designed with ergonomic consideration to maximize donor comfort. Accessories like BP apparatus, blood collection monitor etc.

FEATURES:

- Facility for blood collection from both sides
- Provides a comfortable position for the donor
- Variable positioning for either arm with a Comfortably wide arm-rests
- Tilt adjustment can be done by using remote control
- Comfortable chair type with soft padding for cushioning and rexine cover
- Two Geared Actuator. Ensuring safety and comfort to the donor
- 4 Heavy Duty Lockable castors for easy mobility
- Based on hemodynamic principles. (Donor's head can be lowered immediately and legs can be lifted above heart level so that blood can flow back to the brain and other vital organs, in case of vasovagal attack)

Specification:

Model	BDC-10	Power Consumption	100 W
Power Supply	230 VAC	Length of Arm Rest	60 cm
Weight	60 KG	Width of Arm Rest	15 cm
Lifting Capacity	150 KG	Control	Corded
Size in Normal	620*2050*1440cm(WxDxH)	Remote Upholstery	Soft 4.0 Thick
Size in stretched	620*2050*1100cm(WxDxH)	Movement Actuation	DC motor


 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved