



BER -10 Series Plasma Extractor

| Blood Bank

☰ Details Info

Features:

Plasma Extractor & Plasma Expressor is designed to Extract Blood component from centrifuged Bags. The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag. Its sturdy construction makes it durable and reliable, easy to use and portable. The unit is available in both versions Semi automatic and manual.

Specification:

Model	BER-10	BER-10A
Sensor	-	Infra Red
Clamping Control	Manual	Motor Activated
Pressure Plate	Acrylic	Acrylic