



## ICC-12 ID Card centrifuge

| Blood Bank

☰ Details Info

### ICC-12 ID Card centrifuge

Speed: 800--1800rpm

Centrifugal Force: 54~272g

Capacity: 12 ID-Cards

Digital centrifuge for Gel cards. ICC-12 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time. It can be suitable many types ID-CARDS from different manufacturers.

#### Features

1. Simultaneous centrifuge of up to 12 Gel cards
2. Removable and exchangeable spinning head
3. Pre-defined centrifugal parameters for optimum Gel cards processing: acceleration, speed, braking and time
4. Digital interface visor control and visualization of centrifuge speed and time remaining for processing
5. Programmable audible alarms for end of centrifugal time periods

#### Specification:

Capacity	12 Gel cards
Speed	Two stages (It can be variable according to user need)
Max..Speed	1500rpm
Centrifugal Force	54~272g
Timer	1~20min
Power	AC100-240V, 50/60Hz, 85W
Dimensions	DxWxH=440x340x170 mm
Weight	7.0kg