



ICC-24 ID Card centrifuge

| Blood Bank

☰ Details Info

ICC-24 ID Card centrifuge

Model: ICC-24
 Speed: 500~1500rpm
 Centrifugal Force: 54~272g
 Capacity: 24 ID-Cards

Digital centrifuge for Gel cards. The ICC-24 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time.

It can be suitable many types ID-CARDS from different manufacturers.

Features

1. Simultaneous centrifuge of up to 24 Gel cards
2. Removable and exchangeable spinning head
3. Pre-defined centrifugal parameters for optimum Gel cards processing: acceleration, speed, braking and time
4. Digital interface visor control and visualization of centrifuge speed and time remaining for processing
5. Programmable audible alarms for end of centrifugal time periods

Specification:

Capacity	24 Gel cards
Speed	Two stages (It can be variable according to user need)
Max..Speed	1500rpm
Centrifugal Force	33~294g
Timer	1~20min
Power	AC100-240V, 50/60Hz, 85W