



## ICI-24 ID Card incubator

| Blood Bank

☰ Details Info

### ICI-24 ID Card incubator

Model: ICI-24  
 Temperature: 37C degree  
 Timer: 15min (Variable)  
 Capacity: 24 ID-Cards

Digital incubator for the Gel cards. The ICI-24 ID-incubator holds 2 independent incubation zones with a capacity for 12 Gel cards, and user can set temperature and incubator time. It can be suitable many types ID-CARDS from different manufacturers.

Features :

1. Incubator for 24-Gel Cards in racks at 37°C
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital interface screen for control and visualization of remaining incubation time and real temperature for each zone.
4. Fixed pre-set temperature of 37 degree
5. Pre-set incubation time of 15 minutes, which can be adjust by the user
6. Audible programmable alarm to notify end of incubation period

Capacity	24 Gel cards 2. Temperature: 37°C
Timer	15min (Variable)
Power	AC220V or AC120V 50-60Hz, 300W
Dimensions	DxWxH=380x220x115 mm